**OSSTEM<sup>®</sup>**  
IMPLANT

# 485 KIT

## User Manual

**OSSTEM<sup>®</sup>**  
IMPLANT

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### ■ Description

485 KIT can help safely place a Short Implant on the mandibular molar where the vertical bone mass is insufficient. Its procedure is simple and its blade is rounded in shape to minimize damage to the inferior alveolar nerve. It is a dedicated kit for a Short Implant of TS/US/SS taper fixture, and can help place fixtures of 4 - 8.5mm .

### ■ Indication for use

A device intended for use in dental implant operation.

### ■ Directions for use

1. 485 drill does not have Y-dim of a general twist drill, and the distance to the laser marking refers to the fixture placement depth.
2. There is a 1mm margin above the laser marking to allow for additional drilling depending on the condition of the surrounding bone.
3. The laser marking indicates the drilling depth levels. The bottom line of the markings serves as the baseline.
4. When drilling, move the hand piece perpendicularly in an up and down(pumping) motion.
5. Be sure to spray saline solution during to reduce the friction heat with the bone.

### ■ Contraindications

Patients with the following conditions are contraindicated for surgical procedures

1. Patients with coagulation disorders or bone/wound treatment disturbances.
2. Patients with uncontrolled diabetes, heavy smokers, or patients under the influence of alcohol.
3. Patients with decreased immune function due to chemotherapy or radiation therapy.
4. Patients with intraoral infection or inflammation.
5. Patients with untreatable occlusion/articulation problems and insufficient arch space.
6. Patients ineligible for surgery due to other reasons.

### ■ Caution in use

1. When using Drill Extension, do not apply excessive force.
2. The recommended number of use for the drill and the driver is 50.
3. The maximum permissible torque for 0.9 Hex Driver is 20Ncm and that for 1.2 Hex Driver is 35Ncm.

### ■ Storage and maintenance after use

1. Separate all used tools immediately after the treatment and store them at room temperature after cleaning and drying.
2. Do not leave the instrument in a place vulnerable to contamination.

3. Be sure to sterilize these medical tools in the Autoclave before use.(132°C, 15minutes)
4. Do not use hydrogen peroxide as a disinfectant or cleaner, because it may cause damage or discoloration of TiN Coating, laser marking, color coding, and anodizing.
5. This product comes with a ONE(1) year warranty from the date of its purchase.

### ■ Side Effect

A few problems may occur after the operation (loss of implant stability, damage of prosthesis, etc.). Deficient quality and quantity of the remaining bone, infection, allergic reaction, inferior oral hygiene or uncooperativeness of patient, implant mobility, partial deterioration of tissue, and improper position or arrangement of implants may cause the above mentioned problems.

### ■ General precautions

The surgical technique for a dental implant involves professional and complex procedures. Thus, the relevant professional education and training is required for implant operation.

### ■ Warnings

Inappropriate selection of patients and operation may cause implant failure or loss of supporting bone s. This product should not be used for purposes other than those prescribed by OSSTEM

IMPLANT and should not be altered in any way. Implant mobility, bone loss, and chronic infection may result in the failure of an implant operation.

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DEUTSCHE OSSTEM GmbH