OSSTEM[®]
IMPLANT

Denture 4U KIT

User Manual

OSSTEM[®]
IMPLANT

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15. Multi Angled Abutment 및 Multi Abutment 체결

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■ Product Description

Denture 4U KIT is intended to accurately and stably implant 4 or 6 fixtures per jaw. And it includes guides and drills that help drilling positioning of the drill of up to Ø2.2 for the anterior region and the drill of up to Ø3.0 for the posterior region to implant 4 or 6 fixtures per jaw (based on an edentulous jaw).

■ Intended Use

A device intended for use in dental implant operation.

■ How to Use

※ Before starting the procedure, determine the suitable posterior guide option (0°/17°/30°) on CT and select the appropriate posterior guide. Then, attach it with the anterior guide (the posterior guide can be replaced during the procedure).

1. Anterior Guide 1-point Fixation

- set the position of the jaw center line, and use a mount driver to implant and fix the anchor screw into guide anchor hole of the Denture 4U Guide. guide anchor hole that is located at the exact center of the guide. Ensure to stop the engine when anchor screw touches the guide because the anchor screw could run idle. (※ If the anchor screw is not smoothly implanted, perform drilling first with the anchor drill.)

2. Posterior Guide Positioning

(1st Adjustment of Dental Arch) & Fixation

- After positioning the posterior guide (1st adjustment of dental arch), fix the arch Fixation handle with one hand. Then, use a 1.2 hex driver to fix the arch Fixation screw (1.2 hex).

3. Posterior Guide 1st Adjustment

(Implantation Distance & B/L Angle) & Fixation

- Adjust the implantation distance & B/L angle by sliding and rotating the posterior guide to fit the implantation position. Then, use a 1.2 hex driver to fix the posterior Fixation screw.

4. Ensuring Procedure Safety

- For drilling, there are two ways you can choose to avoid the nervus alveolaris inferior:

- ① Completely open the flap and check the location of the mental foramen. Then, move the indicator to touch the posterior guide's drilling hole, rechecking the location of the mental foramen as well as checking the applicable fixture implantation direction and length.

- ② If you choose not to open the flap, connect the path checker to both sides of the posterior guide's drilling hole and extrapolate the line extending from the path checker based on the panorama or CT to check the location of mental foramen.

5. Posterior Guide 2nd Adjustment

(Implantation Distance & B/L Angle) & Fixation

- After ensuring the procedure safety, perform the 2nd adjustment of the posterior guide. Then, use a 1.2 hex driver to fix the posterior Fixation screw.

6. Posterior Guide 2-point Fixation

- Using a mount driver, implant and fix the anchor screw into the posterior guide anchor hole at each side (left/right). Ensure to stop the engine when anchor screw touches the guide because the anchor screw could run idle. (※ If the anchor screw is not smoothly implanted, perform drilling first with the anchor drill.)

7. Posterior Region Drilling (Ø3.0)

- Drill posterior guide drilling hole with the Ø3.0 drill based on the implantation length. (※ If the distance is short, use the 5-mm drill.)

8. Posterior Guide 2-Fixation Removal

- Using a mount driver, remove the anchor screw on posterior guide.

9. Anterior Guide Positioning

(2nd Adjustment of Dental Arch) & Fixation

- Release the arch Fixation screw (1.2 hex) and arch Fixation handle. Then, position the anterior guide (2nd adjustment of dental arch) and fix again the arch Fixation screw (1.2 hex) and arch Fixation handle.

10. Anterior Region Drilling

(One of 4 Holes Each for Left/Right: Ø2.2 Drilling)

- Drill one of the 4 anterior guide drilling holes for each side with the Ø2.2 drill. (For 6 fixtures per jaw, 2 holes can be chosen among 4 holes located each at the left and right)

11. Denture 4U Guide Removal

- Using a mount driver, remove the anchor screw.

12. Drilling Using Countersink

- Drill with reference to the marked line based on the fixture implantation angle and depth.

- 17° drilling: Line at the bottom of the countersink / 30° drilling: Line at the top of the countersink (※ Using marked line at the top: implant at 1 mm below the bone level / Using marked line at the bottom: implant at the bone level)

13. Drilling Using Taper Drill

- Drill with reference to the sequence appropriate for the applicable fixture that is to be implanted using the taper drill of the Taper KIT.

14. Fixture Implantation

- Using a mount driver, implant the fixture in the posterior/anterior region.

15. Multi-angled Abutment and Multi-abutment Engagement

- Using a 1.2 hex driver, engage the multi-angled abutment with the posterior region and multi abutment with the anterior region.

■ Contraindications

The procedure is contraindicated in patients including:

1. Patients whose blood does not coagulate or who have a problem in bone and wound treatment
2. Patients with uncontrollable diabetes or heavy smokers/alcoholics
3. Patients with depressed immunity due to chemotherapies and radiotherapies
4. Patients with infection or inflammation in the mouth
5. Patients with untreatable occlusion/joint disorder or patients for whom the distance is not sufficient
6. Other patients not suitable for surgery

■ Precautions for Use

1. It is recommended to use the drill and driver 50 times at maximum.
2. Before starting the procedure, it is recommended to establish the procedure plan based on, for example, radiological analysis.

■ How to keep and process after use

1. How to clean and store
 - ① Immediately, separate all used tools after the procedure.
 - ② When visually checked, clean it with running brush in running water until all impurities are removed from the surface. If the structure is complex, wash it several times so that the internal cleaning part is completely cleaned.
 - ③ Rinse the wash under running water.

- ④ after washing, dry at room temperature until completely removed, and store at room temperature.

2. Do not place the product in the place with the risk of infection.
3. Be sure to sterilize these medical tools in the Autoclave After wrapping the KIT with sterile cloth before use. (Temperature: 132°C, Time: 15 min)
4. Do not use hydrogen peroxide that is commonly used in disinfectants and detergents since it can cause damage or discoloration, including TiN coating, laser marking, color coding, and anodizing.
5. To maintain the sterilizing performance of the sterilizer, it is recommended to regularly check the sterilizer through the manufacturer.

■ Side Effects

A few problems (loss of implant stability, damage to the prosthesis, etc.) may occur after the procedure. Such problems may be caused by qualitative and quantitative insufficiency of the remaining bones, infection, allergic reaction, patient's poor oral hygiene or poor cooperation, movement of the implant, local tissue degeneration, or inappropriate location and arrangement of the implant.

■ General Precautions

A surgical procedure for a dental implant is a specialized and complex procedure. Therefore, relevant regular education and training are required for performing an implant procedure.

■ Warning

Inappropriate patient selection and techniques would cause implant failure or loss of supporting bones. This product must not be used for any other purposes than its intended use specified by OSSTEM IMPLANT CO.,Ltd. or modified by any means. Movement of the implant, loss of bones, and chronic infection could cause failure of an implant procedure.

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