Easy Fixture Removal KIT (EFR KIT)

User Manual

OSSTEM[®]
IMPLANT

OSSTEM[®]
IMPLANT

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1. Description

The Easy Fixture Removal Kit consists of Remover Screw, Remover Body, and Torque Wrench for the safe removal of the fixture by reverse rotation.

2. Purpose

The purpose of this product is to remove the fixture without loss of the surrounding bones when the fixture implanted in the alveolar bone needs to be removed due to the failure of its function (fixture damage, connection part damage, severe bone resorption, etc.).

3. Instruction for Use

< Product Guide by Fixture Type >

System	Fixture Size	Product used (Normal/Fracture)	Color
TS	Ø3.5	F3.5	Yellow
	Ø4.0-Ø4.5	F4.0/4.5	Green
	Ø5.0-Ø7.0	F5.0	Blue
SS (compatible with TS)	P4.8	P4.8	Green
	P6.0	P6.0	Blue
US (only for US)	Mini	F3.5	Yellow
	Regular	F4.0/4.5	Green
	Wide	F5.0	Blue

- 1) Check the condition of the fixture to be removed. (*For fixture fracture, refer to the guide below.)
- 2) The TS system has 3 sizes, and the SS system can share products with the TS system; the US system has its own products depending on size.
- 3) Color is divided into yellow, green, and blue depending on 3 sizes.
- 4) The Normal Remover Screw is selected depending on the fixture type and tightened to the fixture using a Screw Driver. **The tightening torque is 80Ncm (However, F3.5 and mini (yellow color) are 60Ncm).** (Check the Torque Wrench scale.)
- 5) As with the previous step, the Remover Body is selected depending on the fixture type and tightened to the Remover Screw counterclockwise by hand. Afterward, turn counterclockwise using a Torque Wrench.
- 6) When the Remover Body comes into contact with the fixture, the rotational torque increases. At this time, the fixture is removed while slowly turning counterclockwise without tipping the Remover Body. (Irrigation on the Remover Body and the fixture is required to prevent heat production.)
- 7) The fixture removed from the bone is fixed to the Fixture Wrench; the Remover Body is then disassembled while being turned clockwise using a Torque Wrench.
- 8) The Remover Screw is released from the fixture while being turned counterclockwise using a Screw Driver.

(User Guide for Fixture Fracture)

- 1) In case of fixture fracture, the Remover Screw can be used as follows:

- 2) Fracture of the TS/SS Fixture can be divided into three cases: ①crack, ②partial fracture, and ③ total fracture of the upper part. 2-1) If crack is formed, or partial fracture is observed, the Normal Remover Screw can be used. 2-2) If total fracture of the upper part is observed (in case of no-connection part), the Remover Screw exclusive for fracture should be used.
- 3) In case of 2-1), the fixture is removed during the removal process in most cases; when only the Remover Screw is released without removing the fixture, however, the Remover Screw should be replaced with the one exclusive for fracture.
- 4) For the US Fixture, the same products can be used for both normal and fracture conditions.
 - * For additional details, refer to catalogues and operation manuals.

4. Precautions

- Familiarize yourself with the instructions before operation.
- This product is a medical device and should be used after sterilization in the autoclave before surgery. (Temperature: 132°C, time: 15 min)
- The Remover Screw and the Remover Body should only be used once. (Re-use of the product results in increased possibility of Remover Screw and Remover Body fracture.)
- This kit cannot be used ①if the Remover Screw is not tightened due to damage to the inner screw of the fixture, ②if there is horizontal fracture in the lower part of the hex or ③if the degree of vertical fracture leaves less than half of fixture; in such case, remove the fixture using a Trephine Drill.
- Dense bone quality and horizontal fixture fracture create difficulty in removal torque. Use Trephine Drill first to remove around 3mm of cortical bone in order to reduce the removal torque before using EFR KIT.
- In case of high torque during fixture removal, irrigate sufficiently to prevent heat production in the Remover Body and the fixture.
- Metal powder from the fixture can be produced due to high torque in the course of fixture removal. Apply sufficient suction during the operation.
- Stop using the system when torque above 400Ncm occurs. Use Trephine Drill to remove around 3mm of cortical bone and resume the use of EFR KIT
- When the remover screw is fractured during fixture removal, remove the fixture using a Trephine Drill.
- May have to use the new Fixture Wrench when separating fixture from Remover Body by Fixture Wrench, due to status of abrasion of the Fixture Wrench.

5. Cautions

- Separate all used tools immediately after surgery. Store them at room temperature after cleaning and drying.
- Do not leave the instruments in a place vulnerable to contamination.
- Do not use hydrogen peroxide as a disinfectant or cleaner (it may cause damage or discoloration of the TiN coating, laser marking and color coding).
- This product comes with a 1-year warranty from the date of purchase.

6. Device information

Refer to the label on the package for device information such as part name, specifications, manufacturing date, lot, and expiration date.

* Note: For general information on implants, see catalogues, operation manuals, or other related materials.

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DEUTSCHE OSSTEM GmbH