Osteo / Osteotome KIT

User Manual

OSSTEM[®]
IMPLANT

OSSTEM[®]
IMPLANT

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■ Product Contents

- 5 Concave-type osteotomes (Φ2.0, Φ2.5, Φ3.0, Φ3.5, Φ4.0)


Designed for maxillary sinus floor elevation for the vertical augmentation of available bone in the maxillary posterior region.

- 5 Expanding-type osteotomes (Φ2.0, Φ2.5, Φ3.0, Φ3.5, Φ4.0)

Designed for the improvement of the initial stability of implants by densifying the trabecula of bone while preserving poor-quality bone instead of cutting.

- Length marking line

Marked in 5mm, 7mm, 8.5mm, 10mm, 11.5mm, 13mm, and 15mm increments. The Laser Marking between 10mm and 11.5mm is indicated as a thick line.

■ How to Use

1) Preparations before use

- Dental implant surgery involves professional and complex procedures. As such, the relevant professional education and training is required for implant surgery.
- This product is a medical device and should be steam-sterilized in the Autoclave before use. (Temperature: 132 °C, time: 15 minutes)

2) How to use

1) Osteo KIT

- Osteo KIT consists of Concave osteotome and Expanding osteotome.
- Concave osteotome: A device used for elevation of the maxillary sinus floor to vertically increase the bone mass of the alveolar bone available in the maxillary posterior region.
- Expanding osteotome: A device used to increase the initial fixation force of the implant, not by bone shaving of the low-quality bone, but by making the bony trabecula dense while preserving the bone.
- Each osteotome has a Stopper to control the depth of the procedure.

2) Osteotome KIT

- Osteotome KIT consists of Concave osteotomes.
- Concave osteotome: A device used for elevation of the maxillary sinus floor in order to vertically increase the bone mass of the alveolar bone available in the maxillary posterior region.
- Each osteotome has a Stopper to control the depth of the procedure.

3) Cleaning and storage after use

- After the procedure, clean the instruments before storage through the following steps.
 1. Immediately separate all instruments used after the procedure and wash them thoroughly with flowing water using a cleaning brush.
 2. Remove moisture completely with a dry cloth or a hot air fan.
 3. After removing the moisture, insert and fix the treatment instruments in the correct position in the KIT case.
 4. After setting is complete, the KIT is wrapped in a surgical bag, sterilized in the autoclave (temperature: 132 °C, time: 15 minutes), and stored at room temperature.

■ Precautions for use

Warnings

As dental implant surgery is highly specialized and involves complex procedures, it is highly recommended to complete a professional training course before the surgery to ensure the safe and effective use of Osstem's implants. Inappropriate patient selection and surgical techniques can lead to implant failure or bone loss. Osstem's implants should only be used for the indicated use. Dental implants should not be changed arbitrarily. Do not use electro-surgical instruments or a laser around Fixtures and Abutments because of the risk of electric shock or burn. Implant mobility, bone resorption, and chronic infection indicate implant failure. If an implant is infected by a patient's bodily fluid for any reason, the implant cannot be used for other patients.

Contraindication

Do not use Osstem's implants if the alveolar bone of the appropriate width and height is absorbed to the extent that it cannot be secured around the implant after implantation. Patients with symptoms such as insufficient bone, poor bone quality, excessive smoking, blood disorder or uncontrolled diabetes may experience lack of osseointegration or failure of the implant.

Precautions for use

Precautions for use and cleaning

1. Osteotomes used after the procedure should be immediately separated from the Stopper, washed, reassembled, sterilized, and stored at room temperature.
2. Do not leave in a potentially infected area.
3. This product is a medical device and should be steam-sterilized in the Autoclave before use. (Temperature: 132 °C, time: 15 minutes)
4. Hydrogen peroxide must not be used as a disinfection and cleaning agent as it may damage or discolor the laser marking.
5. The warranty period for this product is one year after purchase.

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DEUTSCHE OSSTEM GmbH