

Sold More Than 5 Million Units in 70 Countries


OSSTEM[®]
IMPLANT

TS SA Fixture

- Comfortable and fast placement, Stable initial stability
- Various specifications for all indications
- Clinically proven SA surface

Comfortable and Fast Placement Stable Initial Stability

- Excellent self tapping for comfortable placement
- Double threaded structure for fast placement
- Optimized body shape ensures stable initial stability



Various Specifications for All Indications

Straight

TSII (0°)



Taper

TSIII (1.5°)



Taper




TSIV (6°)



Clinically Proven SA Surface

- Crater and micropit formation for bone formation activation
- Optimal surface roughness (Ra 2.5 ~ 3.0 μ m) for stable bone formation
- 97% (3M), 100% (6M) survival rate in early loading cases for maxillary posterior edentulous patients

* Int J Periodontics Restorative Dent.2015;35:21-276



※ Micro Pigs Mandible