



BONE SPREADER KIT (OBSOK)

BONE SPREADER KIT (OBSOK)


- KIT used to expand the ridge of narrowed alveolar bone
- Offset type for convenient procedures
- Components (4 types):
 - OBSO22F, OBSO28F, OBSO35F, OBSO35R



- Use for alveolar bone expansion
- Offset type for easy operation
- Depth marking corresponding to the implant length



Direction for use : refer to the above schematic



Code	Spec	Tip Length*			
		7	8.5	10	11.5
OBSO22F	Thickness	1.15	1.3	1.45	1.6
	Width	2.1	2.2	2.2	2.2
OBSO28F	Thickness	1.15	1.3	1.45	1.6
	Width	2.65	2.8	2.8	2.8
OBSO35F	Thickness	1.3	1.45	1.6	1.8
	Width	3.3	3.5	3.5	3.5
OBSO35R (round type)	Thickness	1.85	2.1	2.3	2.55
	Width	3.3	3.5	3.5	3.5

*unit: mm